



Introducing the CatEyePlug™ for Tissue Culture

The CatEyePlug provides easy transplanting of tissue culture starts and pre-rooted cuttings. The “eye” of the plug allows for the precise and rapid placement of small and delicate plant material, such as tissue culture, micro cuttings and other hard to start plants.

The plug’s unique eye-shaped dibble makes it possible to open the plug with one hand. Once squeezed, the CatEyePlug opens and creates the perfect spot for the placement of your young plant. When you let go of the plug, it closes and holds the tissue culture plant perfectly in position.

Grower Benefits

- Easier and faster planting of tissue culture plants
- Shorter rooting period
- Optimal aeration
- pH balanced
- Even dry down



CatEyePlug™ Technical Data Sheet

DESCRIPTION: The CatEyePlug is utilized for easy transplanting of tissue culture starts and pre-rooted cuttings.

APPLICATION: CatEyePlugs designed for the production of tissue culture plants.

PHYSICAL PROPERTIES:

- PH: 5.7 +/- 0.2
- EC: 0.8 +/- 0.2

FOOD SAFETY: This product is food safe when used in conjunction with your own food safety program.

TO USE:

1. CatEyePlugs arrive pre-moistened. Add water only if needed.
2. Mist seeds/water in cuttings. Plugs should be moist to the touch, not wet.
3. Watering requirements will change based on type of plant and the growing environment. Monitor as needed. CatEyePlugs have excellent water retention characteristics so be sure not to overwater. Before transplanting out of trays, dry plugs down for easy removal

STORAGE & HANDLING INSTRUCTIONS: Store out of direct sunlight in a cool dry location. As with any peat product, wearing gloves when handling is recommended. If product is dry, a dust mask may be worn.

SHELF LIFE: 6 months when stored as directed. Monthly inspections are recommended.

DISPOSAL: Please contact your local waste disposal facility for end-of-life options.